

Our pan-range injection system with a unique audible indicator



Welcome to the Click™!

- 
Patented Haptic Control System
 Increases delivery consistency in C loop* lens profiles
- 
Pre-Loaded
 No handling of the lens
- 
Micro Incision
 2mm tip for micro incision
- 
OVD Free
 Self-lubricating chamber
- 
Multi Product Platform
 Compatible with numerous Lenstec™ products
- 
Double Click™ Indicator
 Provides audible reassurance of successful preparation
- 
Easy Glide Injection
 Compressible injector tip to aid smooth delivery of lens
- 
Slow Release Spring
 Provides a controlled and consistent delivery

Advancing lens delivery in to a new era...

Lenstec™ has taken lens injection beyond the pre-loaded generation, using patented innovation to develop the Click™ system; a pan-range platform to introduce multiple lens designs through the same injection system.

The Click™ Haptic Control System allows delivery of C-Loop* lenses with increased consistency whilst the Easy Glide Tip comfortably condenses when transitioning through the injector, to deliver a smooth consistent motion.

* The patented control system is not necessary for plate haptic design IOLs




Click™ Pan-Range Injection System Guide




STEP 1

Initial click of the double click indicator



STEP 2

Patented haptic control system
Applies to C loop* lenses only



STEP 3

Second click indicates lens folding process is complete



STEP 4

Resulting in smooth delivery of the lens

Completing the injection process

To complete the injection process, the slow-release spring provides controlled and consistent delivery of the lens ensuring its soft tip does not protrude from the cartridge. The Click™ provides zero handling of the implant, with a one-system-fits all platform, reducing the need for theatres to hold multiple products at the same time, eliminating potential errors using the incorrect lens/injector combination. With Click™ it's possible to totally regulate the procedure using the same knife, hand piece, incision, and injector for all compatible lenses.

The Click™ injection system is compatible with the following intraocular lenses...

