



DISPOSABLE OPHTHALMIC SURGICAL KNIVES

Surgeons have spent decades developing different wound construction techniques in order to reduce risk and improve patient outcomes. To this end, choosing the correct, high quality surgical knife is vital. Edgetec™ Disposable Ophthalmic Surgical Knives from Lenstec, offer consistent sharpness and superior quality to enable surgeons to reach these goals and contribute to achieving the best clinical outcomes.



The Product

- Textured for safety and control
- High quality grade 420 stainless steel
- Consistently sharp blade with flawless cutting edge
- Matt finished blade to eliminate reflective glare
- A choice of single or double bevel blade
- Manufactured in the US and fully CE marked
- Available in boxes of 6

The Benefits

- Consistent and uniform incisions
- Superior wound construction
- Self-sealing incisions
- Reduction in healing time
- Reduced risk of infection
- Guarded knives also available



MICRO PHACO SLIT KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE4-1.8	1.8mm	Angled	Full	Up
EDGE4-2.0	2.0mm	Angled	Full	Up
EDGE4-2.2	2.2mm	Angled	Full	Up

PHACO SLIT KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE-2.5A	2.5mm	Angled	Full	Up
EDGE-2.65A	2.65mm	Angled	Full	Up
EDGE-2.75A	2.75mm	Angled	Full	Up
EDGE-2.8A	2.8mm	Angled	Full	Up
EDGE-2.85A	2.85mm	Angled	Full	Up
EDGE-3.0A	3.0mm	Angled	Full	Up
EDGE-3.2A	3.2mm	Angled	Full	Up

SIDEPORT KNIVES MVR DESIGN

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE9-0.8A	0.8mm	Angled	Full	Up
EDGE- 1.00A	1.0mm	Angled	Full	Up
EDGE9-1.2A	1.2mm	Angled	Full	Up
EDGE9-1.4A	1.4mm	Angled	Full	Up

DISC/SPOON KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE10-2.2mm	2.2mm	Angled	Full	Double

STAB KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE-15DEG	15		Full	
EDGE1-30DEG	30		Full	
EDGE1-45DEG	45		Full	

IMPLANT KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE7-4.0	4.0mm	Angled	Full	Up

CRESCENT KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE5-2.0A	2.0mm	Angled	Full	Up
EDGE5-2.0S	2.0mm	Straight	Full	
EDGE5-2.3A	2.3mm	Angled	Full	Up

VITRECTOMY KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE8-20GS	20 GAUGE	Straight	Full	
EDGE8-20GA	20 GAUGE	Angled	Full	
EDGE8-23GS	23 GAUGE	Straight	Full	
EDGE8-23GA	23 GAUGE	Angled	Full	

PRECISION DEPTH KNIVES

MODEL	WIDTH	ANGLE	HANDLE	BEVEL
EDGE0-550	550 Micron		Full	
EDGE0-600	600 Micron		Full	

ALSO AVAILABLE ON REQUEST

CLEAR CORNEAL KNIVES & TRAPEZOID KNIVES

If any sizes are not listed, please contact us.

