




THE PRODUCT

- High quality grade 420 stainless steel
- Consistently sharp with flawless cutting edge
- Matt finished blade to eliminate reflective glare
- A choice of single or double bevel blade
- Available in boxes of 6

THE BENEFITS

- Consistent and uniform incisions
- Superior wound construction
- Self-sealing incisions
- Reduction in healing time
- Reduced risk of infection
- Guarded knives also available

**Stab Knives -**

Allows easy expansion of incisions, able to make a stable sideport and can easily create a port when using a retractor.

**Slit Knives -**

Surgical blade tailored for the initial stages of cataract surgery mainly for use of phaco hand piece.

**Sideport Knives -**

Used to make a secondary incision for the introduction of viscoelastic or aiding in rhexis creation and cortex aspiration and also aids with IOL manipulation.

**Crescent Knives -**

Mainly used for making the scleral tunnel/self-sealing valve incision

**Vitrectomy Knives -**

Offers surgeons exact access points, minimizing the risk of leaks around instruments during posterior segment procedures.

A wider range of knives are available upon request.

Call: +44 1924 382 678 or **Email:** lenstecuk@lenstec.com