

PMMA IOL'S

Technical Specification



LA-501 - Lens Specification

Optic Size:	6.0mm
Optic Type:	Equalconvex
Length:	12.50mm
Haptic Style:	Plate
Angulation:	Step Vault
Positioning Holes:	0
Construction:	1 Piece
Optic Material:	Poly(methyl methacrylate)

A Constant:	115.3
A/C Depth:	3.30mm

Whole Dioptres:	+10.0 to +25.0
Half Dioptres:	+16.0 to +25.0

LS-106 - Lens Specification

Optic Size:	6.0mm
Optic Type:	Equalconvex
Length:	13.0mm
Haptic Style:	Modified C
Angulation:	Step Vault
Positioning Holes:	0
Construction:	1 Piece
Optic Material:	Poly(methyl methacrylate)

A Constant:	118.4
A/C Depth:	5.20mm

Whole Dioptres:	+6.0 to +30.0
Half Dioptres:	+14.0 to +28.0